

X-Cite illuminators have been qualified for the following microscopes:

Microscope Model	EXFO Adapter Part Number
ATTO / BD BIOSCIENCE	
Atto Bioscience	810-00034
Atto Bioscience CARV2	810-00036
CARL ZEISS	
AxioImager	810-00022
AxioObserver	810-00022
Axiophot	810-00022
Axioplan 1, 2	810-00022
Axioskop 2, 2FS, 40	810-00022 (with external device† use 810-00035)
Axiovert 10, 35, 100, 135	810-00035
Axiovert 200	810-00022 (with external device† use 810-00035)
Axiovert 400	810-00035
Axiovert 25, 40	810-00037
Discovery.v8 / .v12 / .v20	No adapter required*
LSM	No adapter required*
Stemi SV11	No adapter required*
Stereo Lumar	810-00022
LEICA	
DM2000, DM6000, DMI4000	810-00033
DMIL	810-00033
DMIRB, DMIRM, DMIRE2	810-00033
DMLB	810-00033
DMLFS	810-00033
DMR, DMRXA, DMRXA2	810-00033
MZFLIII, MZ16F, MZ16FA	810-00038
M10, MZ9.5	810-00038
NIKON	
50i, 80i, 90i	810-00030
AZ100	No adapter required*
Diaphot, Labophot	810-00030*
E600, E800	810-00030
FN-1	810-00030
SMZ1500	810-00030
TE2000, all other "TE"	810-00030
Ti	810-00030
TS100	810-00030*
For "white light TIRF"	810-00045
most other Nikon microscopes	810-00030
OLYMPUS	
BX51, 61	810-00032
CKX41	810-00032
IX51, 71, 81	810-00031
MVX10	810-00046
SZ7, 10, 12, 16	810-00032

† External devices include shutters, filter wheels, etc installed between the adapter and microscope.

* Additional fittings may be required from the microscope manufacturer. Contact EXFO for additional information.

Don't see your microscope model here?

Please contact EXFO to find out if it is compatible with X-Cite! x-cite@exfo.com

April 22, 2008
MG