



The X-Cite 120 has been qualified for the following microscopes:

Microscope Model	EXFO Adapter Part Number
ATTO / BD BIOSCIENCE	
Atto Bioscience	810-00034
Atto Bioscience CARV2	810-00036
CARL ZEISS	
AxioImager	810-00022
AxioObserver	810-00022
Axiophot	810-00022
Axioplan 1, 2	810-00022
Axioskop 2, 40	810-00022 (with external device† use 810-00035)
Axiovert 10, 35, 100, 135	810-00035
Axiovert 200	810-00022 (with external device† use 810-00035)
Axiovert 25, 40	810-00037
Discovery.v8 / .v12	No adapter req'd*
LSM	No adapter req'd*
Stemi SV11	No adapter req'd*
Stereo Lumar	810-00022
LEICA	
DM2000, DM6000, DMI4000	810-00033
DMIL	810-00033
DMIRB, DMIRM	810-00033
DMLB	810-00033
DMLFS	810-00033
DMR, DMRXA, DMRXA2	810-00033
MZFLIII, MZ16F, MZ16FA	810-00038
NIKON	
AZ100	No adapter req'd*
Diaphot	810-00030*
E600, E800	810-00030
FN-1	810-00030
SMZ1500	810-00030
TE2000, all other "TE"	810-00030
TS100	810-00030*
For "white light TIRF"	810-00045
most other Nikon microscopes	810-00030
OLYMPUS	
BX51, 61	810-00032
IX51, 71, 81	810-00031
MVX10	810-00046
SZX12, 16	810-00032

† External devices include shutters, filter wheels and other accessories installed between the adapter and microscope.

* Additional fittings may be required from the microscope manufacturer. Please contact EXFO for additional information.

Don't see your microscope model here? Please contact EXFO to find out if it is compatible with X-Cite! x-cite@exfo.com

June 7, 2007
MG